

Data Sheet 60564	 SteriPack™ <small>CONTRACT MANUFACTURING</small>	Rev 002
------------------	---	---------


General Information	
Trade Name	Sterile Polyester Spun Swab
Common Name	3" Polyester Nasal Swab
Product Code	KXG
Classification	Class 1 (exempt)
Manufacturer	SteriPack USA Ltd LLC 4255 South Pipkin Road, Lakeland, Florida, 33811, USA T: +1 863-648-2333 http://steripackgroup.com/
Contact email	info@steripackgroup.com
Registration #	3009656988
ISO	13485:2016
Cleanroom	ISO Class 8
Product Reference	60564
GS1	860004112716
Quantity per box	2,000
Sterilization	Gamma sterilized Validated at VDmax20 using ANSI/AAMI/ISO TIR13004:2013/(R) 2016, i.e. > 20kGy
Package	
Shelf Life	1 year, from DOM

Swab Information

Manufacturer	US Cotton
---------------------	-----------

Cosmetics

- White polyester bud
- Flexible plastic stick



Functional / Dimensional




- | | | |
|--|---------|-----------------------------------|
| • Bud weight | Minimum | 0.31 grams for 10 Swabs (10 Buds) |
| • Bud diameter | Maximum | 6.1 mm (0.24 in.) |
| • Swab length | Nominal | 73.0 mm (2.87 in.) |
| • Stick diameter | Nominal | 2.5 mm (0.098 in.) |
| • Stick length | Nominal | 70.0 mm (2.756 in.) |
| • Bud fibers keep integrity and adhesion to stick under normal use | | |

Material Type

- 100% Polyester bud
- Polypropylene stick

Package Information

Top Substrate	3.1mil DSR 400 Heat Seal Coated/Ethypel® Performance CP 63gsm
Bottom Substrate	4mil PE Coex Forming Web
Dimension - Width	0.88in / 22.28mm
Dimension - Length	4.33in / 110mm

 CONTRACT MANUFACTURING  SteriPack USA, LTD Lakeland, FL 33811 steripackgroup.com	Sterile Polyester Spun Swab - 60564RevA MFG: EXP: LOT: QTY: 1 <div style="display: flex; align-items: center;"> <div style="border: 1px solid black; padding: 2px; margin-right: 10px;">STERILE R</div> <div style="margin-right: 10px;">⊗</div> <div style="font-size: small;"> (Do not touch bud when removing swab) (No toque la punta al retirar el hisopo) </div> </div>	PULL / TIRAR  PULL / TIRAR
--	--	---

4.33in / 110mm

0.88in / 22.28mm