

Data Sheet 60567	 <b>SteriPack™</b> <small>CONTRACT MANUFACTURING</small>	Rev 001
------------------	---	---------


<b>General Information</b>	
<b>Trade Name</b>	Sterile Polyester Spun Swab
<b>Common Name</b>	3" Polyester Nasal Swab
<b>Product Code</b>	KXG
<b>Classification</b>	Class 1 (exempt)
<b>Manufacturer</b>	SteriPack USA Ltd LLC 4255 South Pipkin Road, Lakeland, Florida, 33811, USA T: +1 863-648-2333 <a href="http://steripackgroup.com/">http://steripackgroup.com/</a>
<b>Contact email</b>	info@steripackgroup.com
<b>Registration #</b>	3009656988
<b>ISO</b>	13485:2016
<b>Cleanroom</b>	ISO Class 8
<b>Product Reference</b>	60567
<b>GS1</b>	860004112723
<b>Quantity per box</b>	2,000
<b>Sterilization</b>	Gamma sterilized Validated at VDmax20 using ANSI/AAMI/ISO TIR13004:2013/(R) 2016, i.e. > 20kGy
<b>Package</b>	
<b>Shelf Life</b>	1 year, from DOM

### Swab Information

<b>Manufacturer</b>	US Cotton
---------------------	-----------

#### Cosmetics

- White polyester bud
- Flexible plastic stick



#### Functional / Dimensional


- Bud weight Minimum 0.20 grams for 10 Swabs (10 Buds)
- Bud diameter Range 3.81 mm - 5.08 mm (0.15in. - 0.20 in.)
- Swab length Nominal 72.0 mm (2.835 in.)
- Stick diameter Nominal 2.6 mm (0.1024 in.)
- Stick length Nominal 70.0 mm (2.756 in.)
- Bud fibers keep integrity and adhesion to stick under normal use

#### Material Type

- 100% Polyester bud
- Polypropylene stick

### Package Information

<b>Top Substrate</b>	3.1mil DSR 400 Heat Seal Coated/Ethypel® Performance CP 63gsm
<b>Bottom Substrate</b>	4mil PE Coex Forming Web
<b>Dimension - Width</b>	0.88in / 22.28mm
<b>Dimension - Length</b>	4.33in / 110mm



0.88in / 22.28mm

4.33in / 110mm