



General Information	
<b>Trade Name</b>	Sterile Polyester Spun Swab
<b>Common Name</b>	3" Polyester Nasal Swab
<b>Product Code</b>	KXG
<b>Classification</b>	Class 1 (exempt)
<b>Manufacturer</b>	SteriPack USA Ltd LLC 4255 South Pipkin Road, Lakeland, Florida, 33811, USA T: +1 863-648-2333 <a href="http://steripackgroup.com/">http://steripackgroup.com/</a>
<b>Contact email</b>	info@steripackgroup.com
<b>Registration #</b>	3009656988
<b>ISO</b>	13485:2016
<b>Cleanroom</b>	ISO Class 8
<b>Product Reference</b>	60577
<b>IFU</b>	An electronic IFU (eIFU) is used for this finished good. See <a href="https://steripackgroup.com/sterile-swabs-manufacturing/">https://steripackgroup.com/sterile-swabs-manufacturing/</a> for instructions for use.
<b>GS1</b>	860004112754
<b>Quantity per box</b>	2,000
<b>Sterilization</b>	Gamma sterilized Validated at VDmax20 using ANSI/AAMI/ISO TIR13004:2013/(R) 2016, i.e. > 20kGy
<b>Package</b>	
<b>Shelf Life</b>	1 year, from DOM

### Swab Information

**Manufacturer** US Cotton



TECHNICAL SPECIFICATION		
<b>STICKS</b>		
MATERIAL	POLYPROPYLENE 9003-07-0	WHITE, HOLLOW
LENGTH	2.7555±.0275 in.	69.3 – 70.7 mm
OUTER DIAMETER	ø.1014±.0069 in.	ø2.40 – 2.75 mm
INNER DIAMETER	ø.0630±.0118 in.	ø1.30 – 1.90 mm
WEIGHT	200±15 gr/1000 Sticks	
<b>FIBER</b>		
MATERIAL	PET, POLYETHYLENE TEREPHTHALATE	POLYESTER
WEIGHT (DRY)	MIN: .20 gr/10 Swabs	
<b>SWABS</b>		
LENGTH	2.8740±.1180 in.	70.0 – 76.0 mm
FIBER HEAD (BUD) DIAMETER	ø.185 ±.035 in.	ø3.80 – 5.60 mm
FIBER HEAD (BUD) LENGTH	.630 +.150/- .132 in.	12.65 – 19.81 mm
FIBER HUMIDITY (INDICATIVE)	N/A	
FIBER ADHESION TO STICK	MIN: 5.00 lbf	
QUANTITY	PRODUCT SPECIFIED	
<b>ADDITIVES</b>		
SOLUTION		
WATER	YES	
BINDER	WE0636	
PRESERVATIVE	NONE	
BACTERICIDE	NONE	
GLUE	N/A	

### Package Information

<b>Top Substrate</b>	3.1mil DSR 400 Heat Seal Coated/Ethypel® Performance CP 63gsm
<b>Bottom Substrate</b>	4mil PE Coex Forming Web
<b>Dimension - Width</b>	0.88in / 22.28mm
<b>Dimension - Length</b>	4.33in / 110mm

<p>SteriPack USA, LTD Lakeland, FL 33811 steripackgroup.com</p>	Sterile Polyester Spun Swab - 60577RevA MFG:                      EXP: LOT: QTY: 1 [STERILE] [R]	
	(Do not touch bud when removing swab) (No toque la punta al retirar el hisopo)	

4.33in / 110mm

Data Sheet 60577



Rev 002

0.88in / 22.28mm

## 1 DOCUMENT HISTORY

Approved By	Status	Date/Time
David Moore	Approved	December 23, 2020 08:34:48 (Eastern Standard Time)
Mirelsa Diaz	Approved	December 23, 2020 10:47:37 (Eastern Standard Time)

Issued By	Status	Date/Time
Mirelsa Diaz	Current	December 23, 2020 16:30:31 (Eastern Standard Time)